

# High-throughput, Cost-Effective Sequencing from BGI and Gencove

The unique combination of DNBSEQ™ sequencing and Gencove's low-pass sequencing software platform enables researchers to efficiently obtain fully-imputed, analysis-ready VCF files from DNA in a single service.

Low-pass whole genome sequencing (LPWGS) is commonly defined as sequencing a genome to an average depth less than 1x coverage. LPWGS, combined imputation, offers a solution to the limitations of genotyping arrays and the high cost of WGS (30x). LPWGS has been shown to provide a substantial increase in statistical power for genome wide association studies (GWAS), and increased accuracy for polygenic risk prediction (PRS) at effective coverages of 0.5x and higher<sup>[1]</sup>.

# LPWGS advantages:

- · Accuracy with less ascertainment bias
- · Complete genomic information to discover novel variation at the population level
- · Fast and high-volume achieved by multiplexing large numbers of samples
- · Easy to set up and future-proofed by updating the reference panel

# LPWGS is the perfect solution for:

- · Biobanks easily and cost-effectively sequence all samples
- · Clinical research ensure accurate results with genome level insight for large cohorts
- · GWAS predict genotype to phenotype associations with greater certainty
- · Pharmacogenomics get fast and expanded profiling

# **Sequencing Service Specification**

BGI's Human Whole Genome Sequencing Services are performed with DNBSEQ™ sequencing technology.



## **Superior Data Quality**

- · PCR and PCR-free library methods are available
- · 100bp/150bp Paired end sequencing
- · Choice of sequencing depth: low pass WGS (1x, 4x, or custom)
- · CAP/CLIA laboratory services available



### **Sequencing Quality Standard**

· Guaranteed ≥80% of bases with quality score of ≥Q30




#### **Turn Around Time**

- · Typically, 18 working days from sample QC acceptance to Data Analysis Report availability
- · Expedited services are available, contact your local BGI specialist for details



# **Project Workflow**

We care for your samples from the start through to the result reporting. Highly experienced laboratory professionals follow strict quality procedures to ensure the integrity of your results.





# DNBSEQ™ Sequencing Technology

DNBSEQ™ is an innovative high-throughput sequencing technology, first developed by BGI's Complete Genomics subsidiary in Silicon Valley. The system is powered by combinatorial Probe-Anchor Synthesis (cPAS), linear isothermal Rolling-Circle Replication and DNA Nanoballs (DNB™) technology, followed by high-resolution digital imaging.

The combination of linear amplification and DNB technology reduces the error rate while enhancing the signal. The size of the DNB is controlled in such a way that only one DNB is bound per active site on the flow cell. This densely patterned array technology provides optimal sequencing accuracy and increases flow cell utilization.



# Sample Requirements


We can process your gDNA, saliva, blood, fresh frozen tissue, cell pellets and FFPE samples, with the following general requirements:

|  | DNA Amount and Concentration                       | Minimum<br>Sample Volume |
|--|--|--------------------------|
| Regular Samples<br>For PCR libraries       | Intact genomic DNA ≥1µg, Concentration ≥12.5ng/µl  | 15µl                     |
| Low Input Samples<br>For PCR Libraries     | Intact genomic DNA ≥200ng, Concentration ≥2.5ng/µl | 15µl                     |
| Samples<br>For True PCR-Free<br>Sequencing | Intact genomic DNA ≥2μg, Concentration ≥12.5ng/μl  | 15µl                     |

# **Gencove's Imputation and Analysis Platform**


## Out-performs genotyping arrays and reduces bias

Gencove's low-pass sequencing at 0.4x and 1x coverage consistently yields higher imputation quality than genotyping arrays in European and African populations, but the difference is especially pronounced in African samples.



Minor allele frequency in sample

### Imputation r^2 in Europeans across technologies



Moreover, based on the experiment design as shown in Table1, low-pass sequencing at an effective coverage of  $\sim$ 0.4x or higher consistently yields more accurate imputed genotype calls at sites of common variation than the Illumina GSA and that this pattern holds across both EUR and AFR cohorts (Figure 1). DNBSEQ<sup>TM</sup> LPWGS performs the best among all the experiments.

Table 1. Details of experiments conducted. Experiments A-D were based on low-pass sequencing, while experiment E used the Illumina GSA v3.0.

| Experiment | Mean Coverage | Samples            | Assay             | Library Prep   | Sequencer       |
|------------|---------------|--------------------|-------------------|----------------|-----------------|
| A          | 0.67          | 120 (3 replicates) | lpSeq             | KAPA HyperPlus | Illumina HiSeqX |
| В          | 1.25          | 120 (3 replicates) | lpSeq             | KAPA HyperPlus | Illumina HiSeqX |
| С          | 1.20          | 1 (30 replicates)  | lpSeq             | KAPA HyperPlus | Illumina HiSeqX |
| D          | 1.26          | 60 (1 replicate)   | lpSeq             | MGlEasy        | BGIseq 500      |
| E(array)   | NA            | 120 (3 replicates) | Illumina GSA v3.0 | NA             | NA              |

Average NRC by non-reference allele frequency for unfiltered sites




Figure 1. Average non-reference concordance for unfiltered SNPs by superpopulation by non-reference allele frequency in 1KGP3.

[1] Li, J. et. al. (2020). Low-pass sequencing increases the power of GWAS and decreases measurement error of polygenic risk scores compared to genotyping arrays. 10.1101/2020.04.29.068452.



# **Request for Information or Quotation**

Contact your BGI account representative for the most affordable rates in the industry and to discuss how we can meet your specific project requirements or for expert advice on experiment design, from sample to bioinformatics.

info@bgi.com www.bgi.com

## **BGI Offices**

# **BGI Americas**

One Broadway, 3rd Floor Cambridge, MA 02142, USA

## **BGI Europe**

Ole Maaløes Vej 3, DK-2200 Copenhagen N, Denmark

## **BGI** Asia

Building NO.7, BGI Park, No.21 Hongan 3rd Street, Yantian District, Shenzhen

## **BGI Australia**

L6, CBCRC, 300 Herston Road, Herston, Brisbane, Queensland 4006, Australia

For Research Use Only. Not for use in diagnostic procedures (except as specifically noted).

Copyright© BGI 2022. All trademarks are the property of BGI, or their respective owners. This material contains information on products which is targeted to a wide range of audiences and could contain product details or information otherwise not accessible or valid in your country. Please be aware that we do not take any responsibility for accessing such information which may not comply with any legal process, regulation, registration or usage in the country of your origin. Note, BGI's genetic testing products have not been cleared or approved by the US FDA and are not available in the USA. For Research Use Only. Unless otherwise informed, all sequencers and sequencing reagents are not available in Germany, USA, Spain, UK, Hong Kong, Sweden, Belgium or Italy. Certain sequencing services are not available in USA and Hong Kong. Please contact a representative for regional availability. The company reserves the right of final interpretation.







We Sequence, You Discover