

**Rapid Identification Kit for B.1.1.7 Lineage (N501Y and P681H) of SARS-CoV-2
(ARMS -PCR)**

For Research Use Only. Not for use in diagnostic procedures.

Instructions for Use*

***Note: Instructions for use are being provided to assist laboratory researchers conducting research related to B1.1.7 lineage (NY501Y and P681h) of SARS-CoV-2. This kit is for research use only, and these instructions should not be used to obtain clinical diagnostic results.**

【Product name】

Rapid Identification Kit for Identifying B.1.1.7 Lineage (N501Y and P681H) of SARS-CoV-2 (ARMS -PCR)

【Package size】 50 tests/kit

【Intended use】

The kit is a qualitative in vitro nucleic acid amplification assay to identify B.1.1.7 lineage (N501Y and P681H) of SARS-CoV-2 from throat swab specimen or sputum confirmed positive SARS-CoV-2 by RT-PCR for research purposes.

Multiple SARS-CoV-2 variants are circulating globally. Several new variants emerged in the fall of 2020, most notably, in the United Kingdom (UK), a new variant strain of SARS-CoV-2 (known as 20B/501Y.V1, VOC 202012/01, or B.1.1.7 lineage) emerged with an unusually large number of mutations. This variant has been detected in numerous countries around the world. This variant has a mutation in the receptor binding domain (RBD) of the spike protein at position 501, where amino acid asparagine (N) has been replaced with tyrosine (Y), N501Y. This variant also has several other mutations, including P681H, near the S1/S2 furin cleavage site, a site with high variability in coronaviruses. Preliminary epidemiologic indicators suggest that this variant is associated with increased transmissibility (i.e., more efficient and rapid transmission).

This kit is intended for research use only and these instructions for use should not be used to obtain clinical diagnostic results.

【Laboratory principles】

The mutation detection is based on allele refractory mutation system (ARMS)-based quantitative PCR (qPCR). Primers and sequence-specific fluorescent probes were designed to tailor key mutation alleles of N501Y and P681H in SARS-CoV-2 genome for wild strain and variant, respectively. ARMS primer can specifically identify and differentiate wild strain from mutation using different fluorescent bands, enabling to differentiate wild strains (501N and 681P) from mutant strains (501Y and 681H) through 2 tests for each specimen in a single run.

In addition, internal reference is included in the kit to monitor the whole reactions, to avoid false negative detection results.

【Key contents】

| Contents(50 tests/kit) | Specification | Quantity | Description |
|--------------------------------|---------------|----------|--|
| Reaction Mix for wild strain | 1 mL/vial | 1 vial | Reagent with primers and probe for amplification of wild strain and internal reference |
| Reaction Mix for mutant strain | 1 mL/vial | 1 vial | Reagent with primers and probe for amplification of mutant strain and internal |

For research use only. Not for use in diagnostic procedures.

| | | | |
|------------------|-------------|--------|--|
| | | | reference. |
| Enzyme Mix | 160 µL/vial | 1 vial | Taq DNA Polymerase, reverse transcriptase, UDG |
| Positive Control | 750 µL/vial | 1 vial | Mix solution of recombinant pseudo-viruses with target genes of wild strain, mutant strain and internal reference. |
| Blank Control | 750 µL/vial | 1 vial | DNase/RNase free H ₂ O |

Notes: Components contained within a kit are intended to be used together. The reagents from a different kit with different lot numbers cannot be mixed.

Materials required but not provided

Table 1 Materials required but not provided

| Item | Validated products |
|--|---|
| Reagent | QIAamp Viral RNA Mini Kit (Cat. No. 52904 for 50 preps, 52906 for 250 preps) by QIAGEN |
| | TIANamp Virus RNA extraction Kit (Cat. No. YDP315-R) by TIANGEN |
| | MGIEasy Nucleic Acid Extraction Kit (Cat. No. 1000020471 for 96 preps, 10000020261 for 1728 preps) by Wuhan MGI Tech Co., Ltd |
| | Solution for sputum decontamination |
| Extraction equipment (Optional) | DNA Sequencing Library Preparation System (MGISP-100RS) by Wuhan MGI Tech Co., Ltd |
| | High-throughput Automated Sample Preparation System (MGISP-960RS) by Wuhan MGI Tech Co., Ltd |
| Consumables | RNase/DNase-free tips for pipettes |
| | Disposable gloves |
| | RNase/DNase-free microcentrifuge tube, 8-tube strips for real-time PCR |

【Storage and shelf-life】

The kit should be stored at temperature lower than -15°C in dark. It is stable with shelf-life for 9 months from date of production in claimed storage condition. Unpacked kit should avoid repeated thaw-freeze (within 6 times). The kit can be transported at temperatures lower than -15°C in dark conditions and is stable for 7 days.

The manufacture date and expiration date are provided on the product labels.

【Applicable instruments】

Applied Biosystems™ QuantStudio 5 Real-Time PCR Systems; Roche LightCycler® 480 Real time PCR System; Fluorescent Quantitative PCR Detection system FQD-96A; Real-Time Quantitative Thermal Cycler MA-6000

【Specimen】

Sample collection

- Collect fresh specimen of throat swabs and sputum. The operation of specimen collection should avoid possible contamination in collection, storage, and transportation. The specimen should be presumed contagious and be operated according to related regulations.
- **Throat swabs:** Carefully take out the swab from package and quickly rotate it around two sides of

For research use only. Not for use in diagnostic procedures.

faucis, throat and tonsil a few times applying pressure to collect as much secretions as possible. Avoid touching tongue. Break the swab stick and put the head into sampling solution in specimen tubes. Screw the tube cap tightly to ensure no leakage.

- **Sputum:** Collect sputum in the early morning after washing mouth. Take a deep breath. Hold the air for a few seconds. Breathe out slowly. Take another deep breath. Cough hard until sputum comes up in mouth. Spit the sputum into the sample bottle. Do this until there is enough sputum to cover the bottom of the bottle. Gas aspiration method can be used to collect sputum for those without sputum. Screw the tube cap tightly to ensure no leakage and seal the tube with film. The sputum should be delivered for testing immediately.

Storage

- The specimen should be kept in proper condition, at temperature lower than -15°C for not longer than 1 week and at temperature lower -70°C for not longer than 6 months.
- Frozen specimen should be thawed thoroughly while avoiding repeated thaw-freeze cycle.

Transportation

The specimen should be shipped in low temperature condition using dry ice or ice bag.

【Laboratory procedures】 (Please read the procedures carefully before your operation)

Sample processing

- The fresh swab specimen should be collected to ensure the qualified RNA in terms of quality and quantity for the assay. RNA should be extracted using Nucleic Acid extracting Kit in line with the manufacturer's instructions. Equivalent volumes of positive control and Blank control should be processed simultaneously. The assay was assessed using the recommended RNA extraction kits by TIANGEN (DP315-R), QIAGEN(50 Preps: 52904, 250 Preps: 52906) and MGIEasy Nucleic Acid Extraction Kit (96 Preps: 1000020471, 1728 Preps: 1000020261) by Wuhan MGI Tech Co., Ltd. 140µL specimen is used by extraction kits from TIANGEN and QIAGEN. 200 µL specimen is needed for kit from MGI to extract nucleic acid manually and more than 160 µL specimen is needed for kit from MGI to extract nucleic acid automatically using High-throughput Automated Sample Preparation System(MGISP-960RS, Cat. No. 900-000154-00) or DNA Sequencing Library Preparation System (MGISP-100RS, Cat. No. 900-000206-00).
- Sputum should be mixed with an equivalent volume of decontamination solution and shaken for about 30 minutes at ambient temperature followed by Nucleic Acid extraction.
- The extracted RNA should be tested immediately or stored at temperature lower than -70°C for test later.

Reagent preparation

- Take out all the kit contents and thaw them thoroughly at ambient temperature. Vortex and centrifuge briefly. The Enzyme Mix should be kept in ice continuously.
- Estimate the number of reactions (N) in the test, which includes the Blank control (1 tube), Positive control (1 tube), and specimens prepared, respectively, for wild and mutant strains.
- Prepare 8-tube strips for PCR based on the estimated N of reactions for wild and mutant strains, respectively. PCR mix 1 for testing wild strain should be prepared as ingredients in table 2 below. Table 3 below was ingredients for preparing PCR mix 2 for testing mutant strain. Pipette 20µL PCR Mix 1 per tube into the 8-tube strips prepared for testing wild strain and 20µL PCR Mix 2 per tube prepared for testing mutant strain, respectively. Cap them tightly and transfer them to sample processing area. The remaining Nucleic Acid Reaction Mix and Enzyme Mix should be stored at temperature lower than -15°C immediately.

Table 2 PCR mix preparation for wild strain of N501Y and P681H

| | Reaction Mix for wild strain (μL) | Enzyme Mix(μL) |
|-----------|-----------------------------------|----------------|
| PCR mix 1 | 18.5×N | 1.5×N |

Table 3 PCR mix preparation for mutant strain of N501Y and P681H

| | Reaction Mix for mutant strain (μL) | Enzyme Mix(μL) |
|-----------|-------------------------------------|----------------|
| PCR mix 2 | 18.5×N | 1.5×N |

Add sample

- Add 10μL extracted RNA from specimens, Blank control, and Positive control respectively into the 8-tube strips prefilled with PCR Mix 1. Cap them tightly and centrifuge them at 2000rpm for 10 seconds.
- Add 10μL extracted RNA from specimens, Blank control, and Positive control respectively into the 8-tube strips prefilled with PCR Mix 2. Cap them tightly and centrifuge them at 2000rpm for 10 seconds.
- Please refer to table 4 below for an example of PCR tube layout in PCR plate.

Table 4 Example of PCR tube in PCR plate

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|----------|
| | PCR-Mix1 | PCR-Mix2 | PCR-Mix1 | PCR-Mix2 | PCR-Mix1 | PCR-Mix2 | PCR-Mix1 | PCR-Mix2 | PCR-Mix1 | PCR-Mix2 | PCR-Mix1 | PCR-Mix2 |
| A | BC | BC | Sample7 | Sample7 | | | | | | | | |
| B | PC | PC | Sample8 | Sample8 | | | | | | | | |
| C | Sample1 | Sample1 | | | | | | | | | | |
| D | Sample2 | Sample2 | | | | | | | | | | |
| E | Sample3 | Sample3 | | | | | | | | | | |
| F | Sample4 | Sample4 | | | | | | | | | | |
| G | Sample5 | Sample5 | | | | | | | | | | |
| H | Sample6 | Sample6 | | | | | | | | | | |

Note: BC-Blank control; PC-Positive control

Real-time PCR

- Set the fluorescent channels: Please refer to the manufacturer’s instructions of thermocycler for detailed information on channel setting.

FAM channel (Reporter: FAM, Quencher: None) for 23063A (501N) in wild strain and 23063T (501Y) in mutant strain of SARS-CoV-2 S-gene

CY5 channel(Reporter: CY5, Quencher: None) for 23604C (681P) in wild strain and 23604A (681H) in mutant strain of SARS-CoV-2 S-gene

VIC/HEX channel (Reporter: VIC/HEX, Quencher: None) for internal reference

Reference Dye: None

Sample Volume: 30μL

- Configure PCR protocol

| Step | Cycles | Temperature | Duration | Fluorescent signal collection |
|------|-----------|-------------|----------|-------------------------------|
| 1 | 1 cycle | 50°C | 5 min | No |
| 2 | 1 cycle | 95°C | 1 min | No |
| 3 | 45 cycles | 95°C | 5 sec | No |
| | | 55°C | 10 sec | Yes |

Data analysis

- Applied Biosystems™ QuantStudio®5 Real time PCR system

Baseline is set as default.

Threshold: In most cases, the auto threshold line function yields satisfying results. In some cases, it can be adjusted manually. In setting threshold, click [Show Plot Setting], select the target gene to view and the “Show: Threshold” as . Adjust the threshold through dragging it by mouse or inputting values directly, then click [Analyze].

- Fluorescent Quantitative PCR Detection system FQD-96A

Baseline is set as default.

Threshold: In most cases, the auto threshold line function yields satisfying results. In some cases, it can be adjusted manually. In setting threshold, click [Analysis settings], select the target gene to view and the “Automatic threshold” as Auto”. Adjust the threshold by inputting values directly, then, click [Save and analyze].

- Real-Time Quantitative Thermal Cycler MA-6000

Baseline is set as default.

Threshold: In most cases, the auto threshold line function yields satisfying results. In some cases, it can be adjusted manually. In setting threshold, click [Analysis], select the target gene to view and the “Automatic threshold” as Auto”. Adjust the threshold by inputting values directly, then, Hit the Enter key on the keyboard.

Quality control

- Blank control: Ct values at FAM and CY5 channels are 0 or no data available. Ct value at VIC/HEX channel is 0, no data available or higher than 38.
- Positive control: Standard curves at channel FAM, CY5 and VIC/HEX channels are all in S-shape with Ct values not higher than 32.
- Above requirements should be met in a single test. Otherwise, the test is invalid and may be retested in line with the package insert.

For research use only. Not for use in diagnostic procedures.

| Quality control metrics | VIC/HEX (Observation) | FAM (Observation) | CY5 (Observation) | Interpretation |
|-------------------------|--|--|--|-------------------------------------|
| Blank Control | No amplification or Ct value is >38. | No amplification | No amplification | Pass. Proceed to sample analysis |
| Positive Control | Sigmoidal amplification curve and Ct value is ≤32. | Sigmoidal amplification curve and Ct value is ≤32. | Sigmoidal amplification curve and Ct value is ≤32. | |
| Blank Control | Sigmoidal amplification curve and Ct value is ≤38. | Sigmoidal amplification curve and Ct value is ≤41. | Sigmoidal amplification curve and Ct value is ≤41. | Failed. Repeat the Run |
| Positive Control | No amplification or Ct value is >32. | No amplification or Ct value is >32. | No amplification or Ct value is >32. | |

【Threshold and reference range】

- Reference range of the kit was determined based on the Receiver Operating characteristic curve and percentile method. Cut-off values for positive 501 and 681 are Ct values lower than 41 in both wild and mutant strains. The identification of wild or mutant strain should be determined by Ct values in combination with ΔCt .
- ΔCt calculated using formula: Ct value in reaction of wild strain- Ct value in reaction of mutant strain.
- Reference range for different allele

| Allele | ΔCt value | Result |
|--------|-------------------|--------|
| 501 | $\Delta Ct > 4$ | 501Y |
| | $\Delta Ct < -4$ | 501N |
| 681 | $\Delta Ct > 4$ | 681H |
| | $\Delta Ct < -4$ | 681P |

b. Cut-off value for internal reference was determined as 38, not higher than 38 as positive.

【Testing result interpretation】

- Standard curves at FAM channels in reactions of both wild stain and mutant strain are in S-shape, ΔCt should be calculated as table 5 below if Ct values are all lower than 41 at FAM channels for both 501Y and 501N or at CY5 channels for both 681H and 681P. In case of no amplification curve in any channel of either wild strain reaction or mutant strain reaction, 45 should be used as Ct value for ΔCt calculation for the reaction without curve. Testing results should be interpreted as table 6 below.

Table 5 ΔCt calculation

| ΔCt | Fluorescent signal | Allele | ΔCt value |
|--------------|--------------------|--------|---------------------------|
| $\Delta Ct1$ | FAM | 501 | $Ct_{wild} - Ct_{mutant}$ |
| $\Delta Ct2$ | CY5 | 681 | $Ct_{wild} - Ct_{mutant}$ |

Table 6 Testing results interpretation by $\Delta Ct1$ and $\Delta Ct2$

| Allele | ΔCt | Result interpretation |
|--------|-------------------|-----------------------|
| 501 | $\Delta Ct1 > 4$ | 501Y |
| | $\Delta Ct1 < -4$ | 501N |
| 681 | $\Delta Ct2 > 4$ | 681H |
| | $\Delta Ct2 < -4$ | 681P |

For research use only. Not for use in diagnostic procedures.

- In case that amplification curve present in only one allele of either wild strain or mutant strain with Ct value not lower than 41, the specimen was in low concentration of virus if Ct value of internal control at VIC/HEX is not higher than 38. The sample with Ct value of internal reference higher than 38 at VIC/HEX should be re-extracted and retested.
- No amplification curve for detecting allele 501 at FAM or 681 at CY5 in both mutation and wild strain was present, the test was valid if Ct value for internal reference at VIC/HEX is not higher than 38. Otherwise, the specimen should be re-tested.

Example for Result interpretation of N501Y and P681H

| Sample | Target (FAM/CY5) | | | | Internal Reference (VIC/HEX) | Interpretation |
|----------|------------------|--|--|----------------------------|---|--|
| | Channel | Wild strain | Mutant strain | ΔCt | | |
| Sample 1 | FAM | Sigmoidal amplification curve and Ct value is <41. | Any | <-4 | Any | 501N |
| | | Any | Sigmoidal amplification curve and Ct value is <41. | >4 | | 501Y |
| | CY5 | Sigmoidal amplification curve and Ct value is <41. | Any | <-4 | | 681P |
| | | Any | Sigmoidal amplification curve and Ct value is <41. | >4 | | 681H |
| Sample 2 | FAM | Sigmoidal amplification curve and Ct value is <41. | Sigmoidal amplification curve and Ct value is <41. | $-4 \leq \Delta Ct \leq 4$ | Any | Possibly contaminated. Reextracted and retest. |
| | CY5 | Sigmoidal amplification curve and Ct value is <41. | Sigmoidal amplification curve and Ct value is <41. | $-4 \leq \Delta Ct \leq 4$ | | |
| Sample 3 | FAM | No Ct or Ct value is ≥ 41 . | No Ct or Ct value is ≥ 41 . | / | Sigmoidal amplification curve and Ct value is ≤ 38 . | Low concentration of virus RNA. |
| | CY5 | No Ct or Ct value is ≥ 41 . | No Ct or Ct value is ≥ 41 . | / | | |
| Sample 4 | FAM | No Ct or Ct value is ≥ 41 . | No Ct or Ct value is ≥ 41 . | / | 0, no Ct or Ct value is >38. | Invalid test. |
| | CY5 | No Ct or Ct value is ≥ 41 . | No Ct or Ct value is ≥ 41 . | / | | Reextracted and retest. |

【Limitation of the assay】

- This kit is intended for research use only for research involving the development of B.1.1.7 lineage (N501Y and P681H strains) of SARS-CoV-2 from positive SARS-CoV-2 throat swab specimen or sputum. These instructions should not be used to obtain clinical diagnostic results.
- Incorrect results can be caused by improper operations in sample collection, transportation or processing, very low concentration of target virus in the specimens, mutations within the viral genome

For research use only. Not for use in diagnostic procedures.

not covered by the kit's primers and/or probe, and uncontrolled external interference factors, such as PCR inhibitor.

【Performance characteristics】

- The package is intact and liquid contents are clear, transparent and no sediments. All contents are in correct quantity as listed in the package insert.
- Positive control is positive at FAM, CY5 and VIC/HEX channel in testing while Blank control is negative at all channels with Ct of internal reference at VIC/HEX higher than 38 or no value.

【Warning and precautions】

- For Research Use Only. Not for use in diagnostic procedures. Please contact BGI Sales in the event of damage to the protective packaging.
- The kit should be stored and transported in claimed conditions. Thaw all kit components thoroughly and centrifuge them briefly before starting an assay. Avoid repeated thaw-freeze cycle.
- All contents in the package are prepared dedicatedly and validated for the intended research purpose. Replacing any of them will affect the performance of the kit. Components contained within a kit are intended to be used together. Do not mix components from different kit lots.
- Separate laboratory areas are recommended for performing predefined procedures of the assay.
 - a) 1st Area: Preparation Area—Prepare testing reagent
 - b) 2nd Area: Sample processing—Process specimen and controls
 - c) 3rd: Amplification Area— Conduct PCR
- All materials used in one area should always be remained in the area and should not be moved or used in other areas. After the assay procedures, the workbench and lab supplies should be cleaned and disinfected timely.
- To prevent the contamination from exogenous RNA, sample addition should follow the sequence of negative control, specimen RNA and positive control. Filtered tips should be prepared and used separately in preparing reagent and sample addition.
- 8-tube strips for real time PCR capped fasten and transferred to specimen processing area immediately after addition of Nucleic Acid Reaction Mix. Ensure to pipette the samples exactly into the reaction mix in PCR tubes and avoid sticking the samples to the inside tube wall. Mineral oil should be added immediately, and the tubes should be capped fasten immediately after the addition.
- After the protocol of amplification is done, remove PCR tubes from the thermal cycler and discard them in a sealable plastic bag for autoclave and decontamination.
- The workbench and lab supplies should be cleaned and disinfected regularly using 75% ethanol or UV light.
- All pipette tips and centrifuge tubes in the assay should be DNase/RNase-free. The used centrifuge tubes and pipette tips should be discarded in waste bin with Clorox (84) disinfectant and disposed with other laboratory wastes after decontamination.
- Operator should receive professional training before operating.

For research use only. Not for use in diagnostic procedures.

【Disclaimer】

- These kits are distributed for research use only, not for clinical diagnostic purposes.

【References】

- [1] LU Rou-jian, ZHANG Ling-lin, TAN Wen-jie, ZHOU Wei-min, WANG Zhong, PENG Kun, RUAN Li. Development and Comparison of Real-Time and Conventional RT-PCR Assay for Detection of Human Coronavirus NL63 and HKU1[J]. CHINESE JOURNAL OF VIROLOGY, 2008(4).
- [2] NIU Pei-hua, LU Rou-jian, LAN Jia-ming, LIU Gao-shan, WANG Wen-lin, TAN Wen-jie. Development of Novel Multiplex Real-time RT-PCR Assays for Detection of MERS-CoV Infection[J]. CHINESE JOURNAL OF VIROLOGY, 2016(3).
- [3] CHEN Yu-jing. Development of two-panel reactions of real-time PCR for detection of 18 types/subtypes of respiratory viruses[D]. 2015
- [4] Starr, Tyler N., Allison J. Greaney, Sarah K. Hilton, Daniel Ellis, Katharine H. D. Crawford, Adam S. Dingens, Mary Jane Navarro, et al. 2020. Deep Mutational Scanning of SARS-CoV-2 Receptor Binding Domain Reveals Constraints on Folding and ACE2 Binding. Cell 182 (5): 1295–1310.e20.
- [5] Gu, Hongjing, Qi Chen, Guan Yang, Lei He, Hang Fan, Yong-Qiang Deng, Yanxiao Wang, et al. 2020. Adaptation of SARS-CoV-2 in BALB/c Mice for Testing Vaccine Efficacy. Science 369 (6511): 1603–7.

【Contact details】



BGI PathoGenesis Pharmaceutical Technology Co.,Ltd.

Building No.7, BGI Park, No.21 Hongan 3rd Street, Yantian District, Shenzhen 518083, China

Manufacturing Site: Jiangsu Macro & Micro-Test Med-Tech Co., Ltd.

Site Address: No.888, Zhujiang Road, Juegang Sub-district, Rudong County, Jiangsu Province

Please Contact: BGI PathoGenesis Pharmaceutical Technology Co.,Ltd.

Customer and Technical Support: info@bgiamericas.com

Website: <https://www.bgi.com/us>